

BS EN 1650:2019

Chemical disinfectants and antiseptics – Quantitative suspension test for the evaluation of fungicidal or yeasticidal activity of chemical disinfectants and antiseptics used in food, industrial, domestic and institutional areas – Test method and requirements (Phase 2, step 1)

Microbiological Solutions Limited (MSL) Gollinrod, Walmersley, Bury, BL9 5NB, UK

Angela Davies, CEO

Customer: Sanglier Ltd Contact name: Hannah Eyres Email: heyres@sanglier.org.uk Address: Lowmoor Buisness Park, Kirby in Ashfield, Nottingham, NG177JZ PO/Quote number: Q002641 Report Date: 28/04/2020 Issue Number: 1

Megan Barrett Laboratory Manager

Peter Thistlethwaite Technical Projects Manager

The test results on this report refer only to the items tested as supplied by the customer. This report shall not be reproduced except in full and with written approval of Microbiological Solutions Ltd. All reports are archived for a minimum of 2 years. The sample will be retained for 1 month unless otherwise requested in writing.

Microbiological Solutions Ltd Gollinrod Walmersley Bury, BL9 5NB



<u>Scope</u>

The standard method BS EN 1650 describes a suspension test for establishing whether a chemical disinfectant or antiseptic has fungicidal/yeasticidal activity. This laboratory test takes into account practical conditions of application of the product including contact time, temperature, test organisms and interfering substance, i.e. conditions which may influence its action in practical situations. The conditions are intended to cover general purposes and to allow reference between laboratories and product types. Each utilization concentration of the chemical disinfectant or antiseptic found by this test corresponds to defined experimental conditions. However, for some applications the recommendations of use of a product may differ and therefore additional test conditions should be used.

Outline of Test Method (Obligatory Test Conditions)

A sample of the test product is diluted in synthetic hard water in products diluted at point of use or distilled water in the case of ready to use products. A test suspension of bacteria and interfering substance is then added to the dilutions and maintained at $20^{\circ}C \pm 1^{\circ}C$ for 15 minutes \pm 10 seconds. At the end of the contact time an aliquot is taken and the fungicidal/ fungistatic activity is immediately neutralised or suppressed by the validated method. The numbers of surviving bacteria in each sample are determined and the reduction is calculated.

The test is performed using *Candida albicans* and *Aspergillus brasiliensis (niger)* as standard organisms.

Acceptance Criteria

The product when tested as above shall demonstrate at least a 4 log₁₀ reduction in viable yeast/fungal counts.

Microbiological Solutions Ltd Gollinrod Walmersley Bury, BL9 5NB



	Deviation	
Name of Product		
Batch Number & Expiry Date	N/S	
Date of Delivery	01/05/2020	
Period of Analysis	06/05/2020-07/05/2020	
Manufacturer / Supplier	Sanglier Ltd	
Storage Conditions	Ambient	
Appearance of the Product	Clear Liquid	
Neutraliser	N6	
Neutralisation Method	Dilution	
Product Diluent	Distilled water	
Test Concentrations	Neat (80%), Mid-range (50%), Non active (0.1%)	
Experimental Conditions	Clean	
Interfering Substance	Clean 0.3g/l Bovine Albumin	
Test Temperature	20°C ± 1°C	
Temperature of Incubation	Fungi – 30°C ± 1°C for 24hr to 48hr	
Identification of the Fungal Strains:	Candida albicans (ATCC 10231)	
Contact Times	5min ± 10s	
Stability and Appearance During Test	No Change Observed	

Deviations from Standard Method

List deviations/There were no deviations from the standard method

Test Result Summary

The test product received has achieved a >4 log reduction against C.albicans when tested under the condition specified in this report.

See page 2 for acceptance criteria and raw data tables below for complete test results.

Microbiological Solutions Ltd Gollinrod Walmersley Bury, BL9 5NB



Validation and controls

Valid	ation suspens	on (Nv₀)	Experim	ental condition	controls (A)	Neutralis	er or Filtratior	Control (B)	Method Validation (C)			
		x =			x =			x =				x =
Vc1	C.a 4	D	Vc1	C.a 45		Vc1	C.a 4	8	Vc1	C.a 42		
		C.a 40			C.a 43			C.a 57			C.a	42
Vc2	C.a 4	D	Vc2	C.a 41		Vc2	C.a 6	5	Vc2	C.a 42		
	$30 \le \overline{x}$ of Nv_0	< 160?		x ¯ of A <u>≥</u> 0.5 Nv0			\overline{x} of B \geq 0.5 M	lv0	x̄ of C≥0.5 Nv0			
	Yes			Yes			Yes			Yes		

<u>Test Results</u>

SOLUTION PROVIDERS				Test C	Conce	entrat	ions (∾/ ∨)						
Test Organism	Suspe	nsion N		Neat			50%			0.1%				
Candida	10^5	188;	161	10^0		0;	0	10^0	0;	0	10^3	158	;	205
albicans	10^6	17;	22		Na	; <	< 2.15		Na ;	< 2.15		Na	;	6.26
ATCC 10231	N_0^{-1}	6.25			R	>	4		R	> 4		٦		-0.01

Microbiological Solutions Ltd Gollinrod Walmersley Bury, BL9 5NB



Test identification Reference: J001683

KEY

No	Log_{10} number of cfu/ml at the beginning of the contact time = $N/10$
Nvo A	is the number of cfu/ml in the validation test suspension at the beginning of the contact time is the verification of experimental conditions control
В	is the neutraliser toxicity control
С	is method validation
Vc	is the colony forming units counted per 1ml of sample
x	is the average of $Vc_1 \& Vc_2$
x wm	is the weighted mean of N
Na	Log ₁₀ number of surviving cfu/ml in the test mixture
R	$(\lg N_O - \lg Na = \lg R)$ is the calculation for reduction in viability
PASS	= lg R greater than or equal to 4
FAIL	= lg R less than 4
>	greater than
≥	equal to or greater than
<	less than
≤	equal to or less than

Microbiological Solutions Ltd Gallinrod Walmersley Bury, BL9 5NB